

Weekly cleaning and inspection record Ultrasound system

Facility Equipment							Location						
						Manufacturer Install date							
Model													
Serial													
Weekly	cleaning and inspecti	on and clear	ning Indica	ate Pass/Fa	il. Should be	performed in	n addition t	to daily safety checks	and infection con	trol.			
		Inspect and clean					Inspect			Test			
Date	QC performed by	Trackball	Monitor	Control panel	Keyboard	Air filters	Power cords	Peripheral cables (if applicable)	ECG cables (if applicable)	Brakes	Wheels	Corrective action? (Y/N)	
Correct	ive action log File cor	rective actio	n service re	eports unde	er the design	ated tab.							
Date	Comments												
Downloa	d more forms at www.lea	rningconnectio	on.com/ultra	sound-accre	ditation-assist								





Weekly cleaning and inspection record Transducer

Facility												
Transduc	er model											
Serial												
Weekly cl	eaning and inspection and clea	uning Indicate Pass/Fail. Sho	ould be performed in additi	on to daily safety checks an	d infection control.							
Date	QC performed by	Cable	Lens	Connector	Case	Corrective action? (Y/N)						
Corrective	e action log File corrective acti	on service reports under the	designated tab.									
Date	Comments											
Download	more forms at www.learningconnect	ion.com/ultrasound-accreditatio	n-assist.									

