



**PHILIPS**

Transducer care

# Maximize your ultrasound investment

Low-cost consultation on the care and  
maintenance of your Philips transducers

Our Philips Clinical Education Specialists value the opportunity to help meet your service needs in this ever-changing healthcare environment. We know it's critical that your ultrasound transducers operate efficiently and that your sonographers are properly trained to today's standards.

With that in mind, we are pleased to offer you a low-cost consultation on the care and maintenance of your Philips transducers, helping to maximize your ultrasound investment.



### Research\* indicates that:

- 4 percent of sonographers do not sheath intracavitary transducers before usage;
- 33 percent of sonographers do not wipe the transducer with cleaning solution after every patient;
- 7 percent of sonographers do not wipe off the gel

\* Kartanginger, R. et. Al. (1997). *Do sonographers practice proper infection control techniques?* Journal of Diagnostic Medical Sonography (13-282-287). Thousand Oaks, California: Sage Publications.

### Transducer care consultation

At the end of this consultation, we will provide a list of activities performed and recommendations or needed actions.

#### This service may help:

- prevent harm to patients;
- reduce costs;
- improve productivity;
- avoid downtime; and
- meet Joint Commission compliance standards.

## For more information

Please call 1 800 522 7022 and speak with our Ultrasound Resource Specialists.