



PHILIPS

Clinical Education

Education News

Philips iE33 Ultrasound System Vision 2012 Basic System Training Curriculum

Our goal at Philips Healthcare is to provide the clinical education you need to make the most of your equipment investment.

With that in mind, we are pleased to offer an e-learning curriculum focusing on the Philips iE33 Ultrasound System. The aim of this comprehensive ten-part series is to provide the Sonographer/Echocardiographer with a comprehensive bundle of courses at an exceptional value.

These ASRT-accredited learning modules are self-paced tutorials consisting of enlightening content and colorful images. There are no prerequisites for this course. You may take all 10 courses or choose the specific ones that suit your needs:

- System Overview
- Annotations and Smart Exam
- 2D Image Optimization
- Color Flow Imaging and Doppler Imaging
- Tissue Doppler Imaging
- Analysis
- Review and Connectivity
- Stress Echo Protocols
- Vascular Imaging
- Contrast Imaging

Purchase the educational series directly from the Philips Learning Center by clicking Register Today. Philips courses are designed to educate customers on the leading edge technologies specific to each system and application. In today's healthcare climate, it is critical to use your ultrasound system to its fullest potential.

Register Today

For questions about classes or purchasing call: 1-800-522-7022



© 2014 Koninklijke Philips N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Please visit www.usa.philips.com/healthcare-medical-education

Printed in The Netherlands.
IMG-14-28888 * Nov 2014