

Fundamentals of Evaluation of Stents with Ultrasound

In today's competitive and dynamic healthcare climate, it is critical to use your medical imaging systems to their fullest potential. Our goal at Philips Healthcare is to provide the clinical education you need to make the most of your equipment investment.



Philips ultrasound university Vascular 202

This Philips Clinical Education full day course for noninvasive vascular diagnosis specialists is developed to meet the climbing need to standardize performance protocols and diagnostic criteria for noninvasive vascular stent evaluations.

The course includes lectures, interactive case study presentations and a hands-on scanning session. Program topics include a brief history of endovascular stent development, national and international research studies, preoperative and postoperative noninvasive vascular testing protocols and diagnostic criteria with case study correlations including DSA angiography, CTA and MRA.

This course includes protocols and diagnostic criteria for abdominal aortic stent grafts, lower and upper extremity arterial stents, peripheral venous stents, carotid artery stents, renal and mesenteric stents, and porto-systemic shunts (TIPS). This course also includes information on the diagnostic criteria used in the diagnosis of stent complications.

Evaluation of Stents with Ultrasound (Vasc 202)



Course objectives.

Upon completion of this course, the learner should be able to:

- Demonstrate knowledge of the anatomy and pathophysiology of different stents into performing the examinations.
- Demonstrate the knowledge of the duplex sonography protocols and techniques for different types of vascular stenting
- Demonstrate the understanding of the anatomy and pathophysiology of different abdominal visceral, peripheral and extracranial cerebral vascular disease.
- Communicate the knowledge of preoperative and postoperative noninvasive vascular testing protocols and diagnostic criteria of different types of arterial occlusive disease and stent grafts.
- Explain normal from abnormal findings and correctly interpret acquired data using B-mode, color and spectral Doppler imaging modalities.
- Increase confidence to incorporate protocols, techniques and interpretation criteria to improve the accuracy of diagnostic and disease management.
- Demonstrate the knowledge of the diagnostic criteria used in the diagnosis of stent complication such as abdominal aortic stent endoleak, stent migration, in-stent re-stenosis and stent fracture.
- Incorporate protocols, techniques and interpretation criteria to improve diagnostic and treatment accuracy.

Prerequisite

Must be a practicing sonographer.

Faculty

Philips clinical education specialists

Locations

Course may be held in Philips central locations in Alpharetta, Georgia; Bothell, Washington; and Cleveland, Ohio. Other locations may be offered.



For more information

Contact a Philips ultrasound clinical services coordinator at 800-522-7022 and visit our education catalog at

www.learningconnection.philips.com/ultrasound

