



PHILIPS

Clinical Education

Smart Exam Protocol Benefits

Our goal at Philips Healthcare is to provide the clinical education you need to make the most of your equipment investment. Virtual instructor-led training (vILT) events use a robust online classroom platform that is specifically designed for highly-interactive, live online learning.

Philips ultrasound general imaging 102VILT

About virtual instructor-led training

Virtual training is a facilitator-led, live online learning event that is delivered in a virtual environment. Participants can be geographically dispersed and also individually connected. Each learner uses their own computer or other compatible device. Virtual training is synchronous, meaning that participants are connected at the same time as the

facilitator. Philips virtual training events typically range from 60-120 minutes in length with a maximum of 10 participants. This socially engaging, purposefully-designed training allows participants the same quality education of an instructor-led classroom without the need or expense of traveling.

Smart Exam Protocol Benefits (GI102VILT)

Course description

This one-hour course will explain the benefits of Smart Exam Protocols.

Course objectives

Upon completion of this course, the learner should be able to:

- Define Smart Exam Protocol
- Describe the benefits of using Smart Exam Protocols
- Understand how to create a Smart Exam Protocol
- Review the SmartExam Quick Guide

Prerequisite

A thorough knowledge and understanding of 2D ultrasound imaging fundamentals and system instrumentation is required for this program.

Audience statement

This course is intended for sonographers, physicians and administrators interested in enhancing workflow, increasing efficiency, improving ergonomics and creating departmental standardizations.

For more information

Contact Philips ultrasound clinical education at 800.522.7022 and visit our education catalog at www.learningconnection.philips.com/ultrasound

