



**PHILIPS**

Clinical Education

## Introduction to MSK Ultrasound of the Wrist and Hand

Our goal at Philips Healthcare is to provide the clinical education you need to make the most of your equipment investment. Virtual instructor-led training (VILT) events use a robust online classroom platform that is specifically designed for highly-interactive, live online learning.

### **Philips ultrasound general imaging 112VILT**

#### **About virtual instructor-led training**

Virtual training is a facilitator-led, live online learning event that is delivered in a virtual environment. Participants can be geographically dispersed and also individually connected. Each learner uses their own computer or other compatible device. Virtual training is synchronous, meaning that participants are connected at the same time as the

facilitator. Philips virtual training events typically range from 60-120 minutes in length with a maximum of 40 participants. This socially engaging, purposefully-designed training allows participants the same quality education of an instructor-led classroom without the need or expense of traveling.

# Introduction to MKS Ultrasound of the Wrist and Hand (GI112VILT)

## Course description

This two-hour Virtual Training course is designed to provide customers basic knowledge to perform musculoskeletal ultrasound of the wrist and hand. Program topics include the indications, technique and sonographic anatomy of the normal wrist and hand.

## Course objectives

Upon completion of this course, the learner should be able to:

- Review anatomy of the wrist and hand
- Discuss indications for a wrist-hand examination
- Demonstrate and discuss patient positioning for wrist-hand evaluation
- Demonstrate and discuss sonographic evaluation of the wrist and hand

## Prerequisite

A thorough knowledge and understanding of 2D ultrasound imaging fundamentals is required for this program.

## Audience statement

This course is intended for physicians and sonographers interested in expanding their knowledge in the application of musculoskeletal ultrasound.

## For more information

Contact Philips ultrasound clinical education at 800.522.7022 and visit our education catalog at [www.learningconnection.philips.com/ultrasound](http://www.learningconnection.philips.com/ultrasound)

